





Burnaid Burn Gel Dressing Pad 10cm x 10cm

\$8.52

In stock
Only **991** left
SKU 807

 Add to Wish List
 Add to Compare



Details

- The Burnaid Burn Gel Dressing consists of a 3mm thick pad of open cell foam impregnated with Burnaid gel. When placed on a burn, a 3mm layer of gel is held on the skin.

From the patient's perspective, they are afforded almost immediate pain relief from the cooling gel and the anaesthetic properties of the Melaleuca Oil. The cooling property of the dressing helps to halt the burn progression and minimizes further tissue damage. Furthermore being in excess of 90% water, the Dressing re-hydrates the burn. Burnaid Dressings are sterile and exhibit mild biocidal properties thus helping to prevent secondary infection.

From a management perspective, Burnaid Dressings provide an effective and simple way of relieving pain and treating burns. The dressing is placed on the burn and stays in-situ for up to two hours. There is no requirement to keep irrigating the burn, nor is there any 'run-off'. Furthermore the need to administer narcotics for pain relief can be greatly reduced. The Dressing doesn't disguise the extent of the burn. When the burn needs to be further examined, the nonadherent dressing simply lifts off and any residual gel rinsed off.

The Burnaid Burn Gel Dressing allows for easy transport to medical attention, as the gel stays 'in-situ' on the skin far better than chilled-water irrigated dressings. Furthermore as the first response treatment in accident and emergency wards or ER's, the Dressing need only be applied once - saving nursing resources and time. The Dressing has many other practical benefits - especially where clean cold water is not readily available.



Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all