



## AEROPLAST Premium Visual Strip 7.5 x 2.5cm Box/100

€19.38

In stock  
SKU AV2000

♥ Add to Wish List

▮ Add to Compare

Extra Wide Strip



### Details

**Size:** 7.5 x 2.5cm **Pricing Unit:** box/100 **Carton QTY:** 36 The AeroPlast Premium Visual plaster is a premium latex free, hypo-allergenic, visual detectable plaster specifically designed for the retail food industry and non-commercial food preparation

### Visually and Metal Detectable

- AeroPlast Premium Detectable has multiple detectable attributes. The first being it's vibrant blue colour, which is the only colour that does not occur naturally in food, making it easily detectable. The second feature is the ability to be detected where metal detection equipment is operational. The ultimate in high security.

### SupaSorb™ Wound Pad

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

### SupaBond™ Adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling. The unique specifications of film and adhesive are key to the products reliability.




ensuring that no matter how wet, steam-filled or hot the environment becomes, this product cannot become separated from the wearer.

### **Premium plaster backing film**

- Aeroplast Visual is made with the ultimate in film technology, providing a bacterial barrier which resists the ingress of water, oil and grease. The standard plaster size is larger for even greater wound isolation.

### **Contraindications / Precautions**

Change dressing daily or more often according to the condition of the wound. As with all adhesive products apply and remove carefully from sensitive or fragile skin. Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician. **This product is not a substitute for direct medical attention from a healthcare practitioner** 



### More Information

|                        |          |
|------------------------|----------|
| <b>MANUFACTURER_MS</b> | Germany  |
| <b>MODEL NUMBER_MS</b> | K61SD521 |
| <b>ECOLABEL_MS</b>     | yes      |