



Non Adhesive Dressing 10 x 10cm

£0.57

In stock
Only **1000** left
SKU 49411

 [Add to Wish List](#)
 [Add to Compare](#)

Details

Melolin consists of a highly absorbent cotton and acrylic fibre pad which is heat bonded on one side to a very thin perforated polyester film. The film side of the dressing is placed next to the wound.

Features: Melolin consists of three layers:

- Low adherent perforated film
- Highly absorbent cotton/acrylic pad
- Hydrophobic backing layer

Benefits:

- Low adherent perforated film
- Allows for rapid drainage of exudate, reducing trauma to healing tissue
- Highly absorbent cotton/acrylic pad
- Patient Comfort
- Cushions and Protects
- Comfortable and minimises pain on removal
- Easy to cut to shape
- Retains its integrity when cut
- Non-sensitising
- Inert and non-sensitising
- Delays time to strike through
- Significant improvement in terms of strike through, therefore requiring fewer dressing changes



Quantity: Single

Dimensions: 10 x 10cm

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all