



Sam Soft Shell Splint 15" Red - 38 x 10.5cm

\$37.19

In stock
Only **1000** left
SKU 490002

 [Add to Wish List](#)
 [Add to Compare](#)

Details

Exceptionally comfortable and versatile, strong and malleable like the Sam Original Splint, the Soft Shell Splint 15" is designed for quick and easy application when a patient needs to wear a forearm splint for more than a few hours or even several weeks.

The soft aluminium core allows it to be moulded to fit the patient's hand, wrist and forearm. With the proper combination of curves, the Sam Splint Soft can be moulded in seconds into any of four basic shapes. The Sam Splint Soft is thicker on one end than on the other, providing two levels of rigidity. Add the comfort of an open-cell foam interior covered in terry cloth and you have the most versatile forearm splint available today.

- Size: 38cm x 10.5cm

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all