



Stethoscope Liberty Classic - Black

\$106.06

In stock
Only **1000** left
SKU 400117

 [Add to Wish List](#)

 [Add to Compare](#)

Details

- Stainless steel chest piece
- Double inner spring binaural
- Non chill ring and diaphragm retaining ring provide excellent patient comfort
- Packaged in a hanging display case

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



Related Products

Check items to add to the cart or select all