



Stethoscope Liberty Classic Tunable Diaphragm - Pink

\$142.73

In stock
Only **1000** left
SKU 400638

 [Add to Wish List](#)

 [Add to Compare](#)

Details

Features:

- Tunable diaphragm
- Innovative design
- Optimum acoustic performance
- Finest craft finish
- Multi patents
- Soft ear tips
- Chest piece with non-chill ring
- Latex Free
- 5 year warranty



- Free name tag with each stethoscope
- Free spare ear tips and diaphragm
- Professionally packaged in a cardboard box

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all