



Thermometer Genius 2 Tympanic

\$984.24

In stock
SKU 400563

 Add to Wish List

 Add to Compare

Details

The Genius 2 thermometer uses its proprietary Peak Select System technology to determine the body's accurate temperature directly from the tympanic membrane.

Features:

- Tapered Probe Design - Facilitates proper alignment with the ear canal.
- Tip Based Sensor Technology - Improves accuracy and repeatability and minimizes the need for recalibration.
- Peak Select System - Records & analyzes up to 100 measurements per scan.
- Rapid Response Time - Easy to Use.
- Programmable Display - Allows the thermometer to be customized according to hospital protocol.
- No-Touch Probe Cover - Enhanced infection control.
- Real Time Temperature - Accurate readings to improve clinical outcomes.
- F°/C° Toggle Button - Customize temperature read outs to follow facility protocol.
- Pulse Timer - Convenient for pulse and respiration monitoring.
- Extended Battery Life - Takes 15,000 readings on 3 AA batteries.
- Easy Load Probe Covers (Thermometer) - Ease of Use.
- Easy Load Probe Covers (Base) - Quick and Easy Loading.
- ASTM and CEN Compliant - Ensures accuracy and dependability.
- FAC System (Fast Ambient Change) - Ensures accurate reading in fluctuating ambient temperatures.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



Related Products

Check items to add to the cart or select all