



Wheelchair Light Weight - Weight 14kgs

\$378.92

Out of stock
SKU 250103

 Add to Wish List

 Add to Compare

Details

Features:

- Lightweight aluminium frame
- Weighs 14kg when ready for transport
- Tipping levers for attendant
- Easy to fold for convenient transport or storage
- Rear wheels 610mm with strong plastic hand rims
- Fixed arm rests and foot rests
- Dual side brakes

Specifications:

- Maximum Load 100kgs
- Seat Width 46cm



- Seat Depth 40cm
- Seat Height 40cm
- Attendant Handle Height 87cm
- Wheelchair Length 100cm
- Wheelchair Width 65cm
- Wheelchair Width Folded 24cm
- Front Wheel 20 cm castor
- Rear Wheel 61cm PU Wheel

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all