



Aqium Hand Sanitiser Ultra

Out of stock
SKU Aqium Hand
Sanitiser Ultra



Product Name

No options of this product are available.

 [Add to Wish List](#)

 [Add to Compare](#)

Hand sanitisers and washing your hands with a soap-free cleanser and water are two of the best ways to help avoid the spread of germs. Aqium Ultra has the added benefit of moisturising your hands with every use.



Details

Hand sanitisers and washing your hands with a soap-free cleanser and water are two of the best ways to help avoid the spread of germs. Aqium Ultra has the added benefit of moisturising your hands with every use.

Benefits

- Aqium Ultra kills 99.99% of germs fast*



- Aqium Ultra contains 70% alcohol
- Aqium Ultra contains Glycerin, Pro-Vitamin B5, Paraffinum Liquidum & Dimethicone to hydrate the skin
- Aqium Ultra provides 25% better skin hydration after first use and up to 40% improvement in skin hydration after multiple uses**
- Aqium Ultra is a quick and convenient way to sanitise your hands on the go
- Aqium Ultra is a hospital strength formulation

* Kills 99.99% of germs according to BS EN 1276:2009.

** Independent skin hydration clinical study (results available on request)

Directions

Apply small amount to palm. Spread and rub over back of hands and finger tips. Allow hands to air dry. Store Below 30 Degrees Celsius. Do not store in direct sunlight.

Ingredients

Alcohol, Aqua (Water), Paraffinum Liquidum, Glycerin, Dimethicone, Panthenol, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Amimomethyl Propanol, Parfum, Butylphenyl Methylpropional, Amyl Cinnamal.