



Aqium Gel Hand Sanitiser 1L

\$16.45

In stock
Only **99** left
SKU 11577

 Add to Wish List

 Add to Compare

Washing your hands properly with water and a soap-free cleanser is one of the best ways to help avoid the spread of germs. But we don't always have access to a cleanser and water. This is where hand sanitisers like Aqium come in handy.



Details

Washing your hands properly with water and a soap-free cleanser is one of the best ways to help avoid the spread of germs. But we don't always have access to a cleanser and water. This is where hand sanitisers like Aqium come in handy.

Benefits:

- Aqium kills 99.99% of germs fast*
- Aqium contains 70% alcohol
- Aqium contains Glycerin and Pro-Vitamin B5 to moisturise the skin
- Aqium is a quick and convenient way to sanitise your hands on the go
- Aqium is non-sticky and air dries quickly
- Aqium is a hospital grade formulation

* Kills 99.99% of germs according to BS EN 1276:2009.



Directions:

Apply 1 - 2 pumps to palm. Spread and rub over back of hands and finger tips. Allow hands to air dry.

Ingredients:

Alcohol, Aqua (Water), Glycerin, Panthenol (Pro-Vitamin B5), Acrylates /C10-30 Alkyl Acrylate Crosspolymer, Aminomethyl Propanol, Parfum, Limonene, Butylphenyl Methylpropional