



## Thermometer Forehead Infrared Pacific

**\$77.00**

Regular Price \$136.73

In stock  
Only **22** left  
**SKU** PACYT1

 Add to Wish List

 Add to Compare

- Fast non-invasive measurement
- Instant one second response time
- Simple one-button operation
- Measurement range 32.0 degrees Celsius- 43.0 degrees Celsius

### Details

## Thermometer Forehead Infrared Pacific

### Features:

- Fast non-invasive measurement
- Instant one second response time
- Simple one-button operation
- Measurement range 32.0 degrees Celsius- 43.0 degrees Celsius
- 0.1 degrees Celsius increment display
- LCD backlight display
- Can be used with or without sound
- Auto shut-off after 30 seconds
- Requires 2 x AA batteries
- Accuracy:
  - 32.0 degrees Celsius- 34.9 degrees Celsius: +/- 0.3 degrees Celsius.
  - 35.0 degrees Celsius- 42.0 degrees Celsius: +/- 0.2 degrees Celsius.
  - 42.1 degrees Celsius- 42.9 degrees Celsius: +/- 0.3 degrees Celsius.



## **Important:**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare providers. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare providers before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



## More Information

**SPECIAL PRICE**

\$70.00