



AEROCREPE Light Cotton Crepe Bandage 5cm x 4M Wrap/12

\$35.96

In stock
SKU ACL50

 Add to Wish List

 Add to Compare

Light Crepe Bandage 5cm




Details

Size: 5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 40 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

Light Crepe

AEROCREPE light cotton crepe bandages are of a standard conforming nature which is durable and suitable for dressing retention and light support requirements.

Easy to Open Wrapper

AEROCREPE light cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes