



## AEROCREPE Medium Cotton Crepe Bandage 2.5cm x 4M Wrap/12

**\$25.49**

In stock  
**SKU** ACM25

 Add to Wish List

 Add to Compare

Medium Crepe Bandage 2.5cm



### Details

**Size:** 2.5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 60 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

### Medium Crepe

AEROCREPE medium cotton crepe bandages are of a standard conforming nature which is durable and suitable for dressing retention and medium support requirements.

### Easy to Open Wrapper

AEROCREPE medium cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



### More Information

<b>MANUFACTURER</b>	Germany
<b>MODEL NUMBER</b>	K61SD521
<b>ECOLABEL</b>	yes