



Mobile Base for Magnifying Lamp

\$1,621.20

In stock
SKU ACR05-MB

 Add to Wish List

 Add to Compare

Mobile Base for Magnifying Lamp



Details

Size: 93 x 50 x 50cm **Pricing Unit:** unit/1 **Carton QTY:** 1

Magnifying Lamp

A magnifying lamp is integral to accurately diagnose and treat small injuries and provide sufficient light for casualty treatment. The mobile base fits onto the Magnifying Lamp ([ACR05-ML](#)).



More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes