



## AEROINSTRUMENTS Stainless Steel Tweezers 7.5cm

**\$2.99**

In stock  
**SKU** AF9



 Add to Wish List

 Add to Compare

Tweezers Stainless Steel



### Details

**Size:** 7.5cm **Pricing Unit:** unit/1 **Carton QTY:** 2000 The AEROINSTRUMENT range offers a wide range of surgical instruments manufactured from high quality stainless steel components to suit first aid requirements. Assorted sizes and specifications have been selected for use in common first aid procedures.



## More Information

<b>MANUFACTURER</b>	Germany
<b>MODEL NUMBER</b>	K61SD521
<b>ECOLABEL</b>	yes