



AEROFIX Zinc Oxide Tape 2.5cm x 5M Box/12

\$25.72

In stock
SKU AFT25



 Add to Wish List

 Add to Compare

Zinc Oxide Tape on spool with cap




Details

Size: 2.5cm x 5M **Pricing Unit:** box/12 **Carton QTY:** 30 AeroFix Zinc Oxide Tape is suitable for securing dressings to the skin and strapping joints and sporting injuries. **FEATURES & BENEFITS**

Individual Spool

AeroFix Zinc Oxide Tape comes in an individual spool with cap ensuring the tape is kept in optimum condition.

Superior Strength

AeroFix Zinc Oxide Tape is well suited more heavy duty applications such as strapping sporting injuries, strains and sprains. 



More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes