




AEROCASE Large Metal Cabinet 41 x 61 x 18cm

\$194.48


In stock
SKU AKM004

 Add to Wish List

 Add to Compare

Metal Cabinet Large Side Opening

Details

Size: 61 x 41 x 18cm **Pricing Unit:** cabinet **Carton QTY:** 1 The AeroCase range of wall mountable cabinets cater to a broad range of first aid kit requirements in the industry, and includes both metal and plastic constructions with multiple sizes and configurations. Cabinet size should reflect the size and capacity requirement for the specific workplace or industrial kit. 



More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes