



AEROCASE Red and Black Plastic Tacklebox with 2 Trays 20 x 40 x 23cm

\$69.77

In stock
SKU AKP002BR

 [Add to Wish List](#)

 [Add to Compare](#)

Plastic Case Red and Black 2 Tray



Details

Size: 20 x 40 x 23cm **Pricing Unit:** case **Carton QTY:** 6 The AeroCase Range of plastic cases caters to a broad range of first aid kit requirements in the industry.

Toolbox Style Plastic Cases

These top quality plastic toolboxes are made in Australia and feature lift-out or cantilever trays so that all products are visible and readily accessible.



More Information

| | |
|---------------------|----------|
| MANUFACTURER | Germany |
| MODEL NUMBER | K61SD521 |
| ECOLABEL | yes |