



## AEROCASE Large White Plastic Cabinet with Key Latch 32 x 37 x 18cm

**\$94.12**


In stock  
**SKU** AKP003

 [Add to Wish List](#)

 [Add to Compare](#)  
Plastic Cabinet White



### Details

**Size:** 32 x 37 x 18cm **Pricing Unit:** cabinet **Carton QTY:** 1 The AeroCase range of wall mountable cabinets cater to a broad range of first aid kit requirements in the industry, and includes both metal and plastic constructions with multiple sizes and configurations. Cabinet size should reflect the size and capacity requirement for the specific workplace or industrial kit. This cabinet has a key latch. 



### More Information

|                     |          |
|---------------------|----------|
| <b>MANUFACTURER</b> | Germany  |
| <b>MODEL NUMBER</b> | K61SD521 |
| <b>ECOLABEL</b>     | yes      |