



## AEROSWAB Sterile White Non-Woven Swab 10 x 10cm (Packs of 3)

**\$11.78**

In stock  
**SKU** ANS103S

 [Add to Wish List](#)

 [Add to Compare](#)

Non-Woven Swab - 3pcs



### Details

**Size:** 10 x 10cm **Pricing Unit:** box/25 **Carton QTY:** 24 The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

### Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.

### Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application.

### Sterile

- AEROSWAB™ gauze swabs are also available in a sterile form for minimising infection risk when dressing clean wounds





### More Information

<b>MANUFACTURER</b>	Germany
<b>MODEL NUMBER</b>	K61SD521
<b>ECOLABEL</b>	yes