



AEROTAPE White Microporous Paper Tape 2.5cm x 5M Box/12

\$14.27

In stock
SKU AP325



 [Add to Wish List](#)

 [Add to Compare](#)

[Microporous Paper Tape](#)



Details

Size: 2.5cm x 5M **Pricing Unit:** box/12 **Carton QTY:** 30 AEROPORE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroPore Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.

Hypoallergenic

- AEROPORE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Individually Wrapped

- AEROPORE Microporous Paper Tape is available in boxes of individually wrapped rolls ensuring the tape is kept in optimum condition.





More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes