



## AEROTAPE Tan Microporous Paper Tape 2.5cm x 5M Box/12



**\$14.27**

In stock  
**SKU** AP325T

 Add to Wish List

 Add to Compare

AEROTAPE Microporous Tan Tape



### Details

AEROTAPE Microporous Tan Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroTape Microporous Tan Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.

### Hypoallergenic

- AEROTAPE Microporous Tan Tape has a very low irritant risk and is suitable for use on fragile and delicate skin.



### More Information

|                     |          |
|---------------------|----------|
| <b>MANUFACTURER</b> | Germany  |
| <b>MODEL NUMBER</b> | K61SD521 |
| <b>ECOLABEL</b>     | yes      |