



AEROTAPE White Microporous Paper Tape 5cm x 5M Box/6

\$19.61

In stock
SKU AP350

 Add to Wish List

 Add to Compare

Microporous Paper Tape



Details

Size: 5cm x 5M **Pricing Unit:** box/6 **Carton QTY:** 30 AEROPORE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroPore Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.

Hypoallergenic

- AEROPORE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

Ease of Use

- AEROPORE Microporous Paper Tape is available in a ready to use dispenser for easy dispensing and tear off

Individually Wrapped

- AEROPORE Microporous Paper Tape is available in boxes of individually wrapped rolls ensuring the tape is kept in optimum condition.





More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes