



AEROPAD Non-Adherent Dressing 10 x 20cm Box/50

\$52.77

In stock
SKU APD102S

 Add to Wish List

 Add to Compare

Low Adherent Dressing Pad Box/50



Details

Size: 10 x 20cm **Pricing Unit:** box/50 **Carton QTY:** 8 The AEROPAD Non-Adherent dressing is recommended for treating cuts, grazes and minor burns using specialized SupaSorb technology.

Premium Absorbency

- Premium absorbent inner material for weeping and bleeding wounds

Non-adherent

- Special Non-adherent film on each side of the pad minimizing wound trauma when changing dressings.

Individually packed and sterile

- Dressings are available individually on demand and sterile for wound cleanliness.



More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes