



## 25G Needle, 38mm Box of 100

**\$15.33**

In stock  
**SKU** ASN25G38

 [Add to Wish List](#)

 [Add to Compare](#)

25G Needle, 38mm Box of 100

### Details

Useful for vaccination of obese patients, the 25G needle is sterile and has great flow rate.



## More Information

|                     |          |
|---------------------|----------|
| <b>MANUFACTURER</b> | Germany  |
| <b>MODEL NUMBER</b> | K61SD521 |
| <b>ECOLABEL</b>     | yes      |