



AEROPROBE Double-Ended Splinter Probe 11cm

\$2.00

In stock
SKU ASP200



 Add to Wish List

 Add to Compare

Reusable Splinter Probe Double-Ended



Details

Size: 11cm **Pricing Unit:** unit/1 **Carton QTY:** 5000 Double ended splinter probe for removal of deep splinters.



More Information

| | |
|---------------------|----------|
| MANUFACTURER | Germany |
| MODEL NUMBER | K61SD521 |
| ECOLABEL | yes |