




## Disposable Clear Plastic Kidney Dish 700mL (230mm)

**\$2.02**

In stock  
**SKU** KD230

 Add to Wish List

 Add to Compare  
Plastic Kidney Dish

### Details

The AEROSUPPLIES stationary range covers a number of common stationary requirements applicable to the First Aid Industry.

#### Disposable Clear Kidney Dish

Kidney dishes are a widely used product in surgical procedures for holding small instruments and more importantly for rinsing with sterile water. The dishes are also used to contain samples for laboratory analysis.



### More Information

<b>MANUFACTURER</b>	Germany
<b>MODEL NUMBER</b>	K61SD521
<b>ECOLABEL</b>	yes