



## HEARTSINE Samaritan 350P Semi-Automatic Defibrillator

**\$2,367.20**

In stock  
**SKU** PAD-350P

 Add to Wish List


 Add to Compare

Defibrillator Samaritan PAD 350P

### Details

**Size:** - **Pricing Unit:** unit/1 **Carton QTY:** 8

- 8-year warranty
- World Leading CPR Advisor
- Semi-Automatic Shock Delivery
- Easy to Follow Visual and Verbal Guides
- Durable
- Advanced Technology
- Single Battery/Electrode Cartridge

 The HeartSine Samaritan 350P **Defibrillator** has all the features of the [500P Defibrillator](#) except for the CPR advisor technology. This public access **defibrillator** has a 8-year warranty and a unique combined battery and pad-pack with a 4-year life from date of manufacture. After an event the complete pad-pak is replaced, which means you have a brand new battery for the next event. This is more cost-effective than many other devices with expensive so-called long-life batteries when costed over the life of the device. For use in hotels, [schools](#), [restaurants](#), companies, clubs, shopping centres, sports clubs, aviation and shipping, doctor's offices, hospital. The PAD 350P **Defibrillator** device can be used for SCA victims from the age of 8 and above utilising the Adult Pad Pak. For smaller children, the **PAD 350P** has a Special Paediatric Pad Pak available. [HeartSine](#) offers the only device to incorporate the Electrodes and Battery into one cassette for ease of use and the reduction of downtime following deployment.



- 8-year warranty - The longest in the industry
- Free replacement Pad-pak after genuine use (conditions apply)
- IP56 Dust and Moisture ingress protection - the highest in the industry
- Metronome for correct speed CPR
- Software updates via USB
- Made in Ireland by the inventors of the first portable **defibrillator**

The PAD 350P Defibrillator is maintenance-free as are all HST devices. The device is fully CE Marked and FDA approved for use around the world. Software updates are free of charge. The hardware is guaranteed for 8 years and the unit is maintenance-free with only one replaceable part after use. The HST samaritan PAD model 350P: suitable for use by trained lay users, semi-professionals and professional users alike. The device is a lifesaver for private, business, commercial and professional use. [Register your Product Warranty Here.](#)



### More Information

<b>MANUFACTURER</b>	Germany
<b>MODEL NUMBER</b>	K61SD521
<b>ECOLABEL</b>	yes