



AERORESCUE Collapsible Basket Aviation Stretcher

\$1,982.20

In stock
SKU STR-09



 Add to Wish List

 Add to Compare

Basket Aviation Stretcher Collapsible



Details

Size: 220 x 64.5 x 16.5cm **Pricing Unit:** unit/1 **Carton QTY:** 1 AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-09

Weight limit: 205kg The collapsible basket stretcher has the same application as the STR-08 but detaches in half for ease of storage. Folds to 126cm x 64.5cm x 21cm



More Information

| | |
|---------------------|----------|
| MANUFACTURER | Germany |
| MODEL NUMBER | K61SD521 |
| ECOLABEL | yes |