





AEROINSTRUMENTS Stainless Steel Universal Shears with Plastic Tip 15cm

\$2.96

In stock
SKU USP15-P

 Add to Wish List
 Add to Compare
Universal Scissors



Details

Size: 15cm **Pricing Unit:** unit/1 **Carton QTY:** 480 The AEROINSTRUMENT range offers a wide range of surgical instruments manufactured from high-quality stainless steel components to suit first aid requirements. Assorted sizes and specifications have been selected for use in common first aid procedures.



More Information

MANUFACTURER	Germany
MODEL NUMBER	K61SD521
ECOLABEL	yes