



## SebiRinse (Conditioner) - 200g

**\$12.05**

In stock  
Only **96** left  
**SKU 11783**

 Add to Wish List

 Add to Compare

Protective Conditioner



### Details

Specially formulated for use after scalp cleansing and anti-dandruff treatments, this gentle conditioner protects your hair and helps prevent dryness. SebiRinse restores shine and adds body, leaving your hair feeling great.

### Benefits:

- Visibly rejuvenates and protects treated hair.
- Removes any lingering odours from hair that's been treated with therapeutic shampoos.
- Non-oily, water-based formulation avoids waxy build-up on hair, which can aggravate scalp conditions.
- Will not damage permed or coloured hair.
- Soap-free and gentle enough to use every day.

### Directions:

After using Sebitar, massage Sebirinse through hair. Leave on for approximately 2 minutes. Rinse.



## Ingredients:

Aqua (Water), Cetearyl Alcohol, Glycerin, Behentrimonium Methosulphate, Phenethyl Alcohol, Panthenol, Ceteth-20, Parfum (Fragrance)