



## La Grande Large Spacer

**\$24.15**

In stock  
Only **50** left  
**SKU** BHC-LG-1

 Add to Wish List

 Add to Compare



The e-chamber La Grande is a spacer for effective delivery of medication. The high capacity chamber (510ml) is made from durable polypropylene plastic, which is entirely recyclable. This is also anti-static to ensure minimum medication adheres to the chamber, unlike alternatives. A perfect option for preventers to ensure maximum drug delivery.

### Details

Large format asthma spacer ideal for older children and adults.

Features include:

- 510ml
- Anti-static plastic with embedded particles - can't be scratched or washed off
- Single patient re-usable
- Easy cleaning
- Universal inhaler fitting
- 15mm-22mm round mouthpiece
- 100% recyclable
- Inhaler can be conveniently stored inside
- Retail packed and available in quantities of 1, 10, 25 or 50 per box