




## FORTRESS 3-IN-1 COMBO TENS, EMS & MASSAGE

**\$214.50**

In stock  
Only **1000** left  
**SKU** FORTRESS  
2420

 [Add to Wish List](#)

 [Add to Compare](#)

The Fortress 3-in-1 is a combination TENS, EMS & Massage device.

- Classic TENS incorporates a wide variety of output combinations ensuring more effective treatment
- Five modes: B(Burst), N(Normal), M(Pulse Rate Modulation), SD1(Pulse Duration Modulation) & SD2(Pulse Duration Modulation)



### Details

Classic TENS incorporates a wide variety of output combinations ensuring more effective treatment

Five modes: B(Burst), N(Normal), M(Pulse Rate Modulation), SD1(Pulse Duration Modulation) and SD2(Pulse Duration Modulation)

Timer of 15 ,30, 60 minutes (or 15,30,Continuous)

It can be used on a wide range of chronic and acute pain conditions: Including cervical, lower back, lumbago ,leg pain, sciatica, arthritis, tendonitis, fractures, tennis elbow and muscle and joint pain