



ZOLL AED 3 Trainer

\$803.00

In stock
SKU
8028-000001-07

 [Add to Wish List](#)

 [Add to Compare](#)

A training unit that closely mimics the clinical ZOLL AED 3 so students get a realistic rescue experience, with Real CPR Help feedback. For training only.

Details

The ZOLL AED 3 Trainer lets CPR and first-aid instructors give students a realistic AED experience without a clinical device. It closely mirrors the look and prompts of the clinical ZOLL AED 3, and includes Real CPR Help so trainees get real-time compression feedback during practice.

Best for: training organisations, workplace first-aid trainers, schools and clubs running CPR sessions.

Includes:

- CPR Uni-padz electrode training harness
- Trainer electrode gels (one pair) and trainer electrode liner (one pair)
- Batteries (4 x D-Cell)
- Trainer operator's guide
- 12-month factory warranty

Features: Real CPR Help, pads-on-manikin detection, active LCD display, on-trainer configuration and USB software updates.



Important: this is a training device only. It does not analyse a real heart rhythm and cannot deliver a shock. It must never be used on a person in an emergency.