



ZOLL AED Pro ECG Electrodes (box of 600)

\$333.30

In stock
SKU 8900-0003

 [Add to Wish List](#)

 [Add to Compare](#)

Genuine ECG monitoring electrodes for the ZOLL AED Pro, supplied in a box of 600.



Details

Genuine ZOLL ECG electrodes for the AED Pro's 3-lead ECG monitoring function. Supplied in bulk (600 per box) for training organisations and clinical sites that monitor regularly.

Best for: AED Pro sites using the device's monitoring capability.

Includes 600 ECG monitoring electrodes. Compatible with the ZOLL AED Pro monitoring lead set.