



WM1 Workplace First Aid Kit Wallmount - Plastic Case by Trafalgar

\$240.34

In stock
Only **998** left
SKU 876479

 Add to Wish List

 Add to Compare

This case is made from sturdy ABS plastic and contains an 'O' ring to help resist dust and water.

The case includes shelving and compartments for easy and hassle-free storage.

Covers your workplace safety needs.



Details

Wall Mount Workplace First Aid Kit: Plastic Case

Easy access to your first aid equipment is essential in your workplace!

Keep your workplace safe and covered for medical emergencies by having a mounted first aid kit in an easily accessible area.

This case is made from sturdy ABS plastic and contains an 'O' ring to help resist dust and water.

Products are stored in shelving and compartments for easy and hassle-free storage.

Covers your workplace safety needs.



Dimensions:

300 x 390 x 150mm (W x H x D)

Content

Burns, Bites, Stings

- Antiseptic Spray 50ml x1
- Hydrogel Sachet 3.5g x5

Eye Injuries

- Eye Pad x4
- Sodium Chloride (saline) Pods 20ml x8

CPR, Rescue, Survival

- Blanket - Thermal x1
- Resuscitation-Safe Face Shield - Disposable x1

Essentials

- Adhesive Non-Woven Tape 2.5cm x 5cm x1
- Biohazard Bag x1
- First Aid Booklet x1
- Forceps - Disposable x1
- Nitrile Gloves Pack (pairs) x5
- Pen x1
- Plastic Bag Zip Lock (small) x2
- Plastic Bag Zip Lock (medium) x2
- Plastic Bag Zip Lock (large) x2
- Register of Injuries Book x1
- Safety Pins - Assorted x12
- Scissors x1
- Splinter Probe - Disposable x10