



## 3M 2330 Red Dot Resting EKG Electrodes (100)

**\$33.58**

In stock  
Only **1000** left  
**SKU** 400261

 Add to Wish List  
 Add to Compare



### Details

#### Description

The 3M 2330 Red Dot Resting EKG Electrodes are designed to make removal painless and easy, yet remain sticky even after multiple repositionings made. It comes with a tab that allows the lead wire placement to be much more convenient.

It ensures close and intimate skin contact to help in accurate EKG traces, even on challenging sites. Plus, the semi-stiff back prevents the electrode from rolling and lifting off.

#### Specifications

Suitable for Critical Care, Emergency Medical Service (EMS), Emergency Room (ER), Intensive Care Unit (ICU), Pediatric, Physician Office, and Operating Room (OR).

Electrode Size: 1.2" x 0.8"; 2cm x 2cm

Type: Resting

Gel Type: Conductive Adhesive

Radiolucency: Radiolucent

Skin Type: Fragile to Normal

Material: Plastic Film

Quantity: 100pcs/pack



---

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all