



## Active Seal Pediatric Arm Cast Bandage Protection

**\$45.89**

In stock  
Only **1000** left  
**SKU** 20363

 [Add to Wish List](#)

 [Add to Compare](#)

### Details

Active Seal Pediatric Arm Cast Bandage Protector can be used for swimming, water sports, whirlpools, hot tubs, showering or bathing to protect casts and bandages.

- It is a safe, economical protector for casts and bandages during recreational water activities or daily showering and bathing
- Easy application eliminates the need for pumps, hook & loop closures and adhesive tapes that may irritate skin
- Made of TPE (a high performance thermoplastic elastomer that stretches easily)
- Maintains a waterproof seal
- Reusable and folds compactly
- With normal care, lasts for three to four weeks of daily usage. (May be affected by chemicals and heat)
- Length 48cm

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.