



AlcoCheck FC300E Fuel Cell Hand Held Breathalyser

\$632.50

In stock
Only **1000** left
SKU 400590

 [Add to Wish List](#)

 [Add to Compare](#)

Alcohol Testing: AlcoCheck™ FC300e Fuel Cell hand held breathalyser.

Details

- Alcohol Testing: AlcoCheck™ FC300e Fuel Cell hand held breathalyser.
 - Evidential breath alcohol tester
 - Cumulative test counter
 - Onboard memory of 1000 tests
 - Portable printer capable
 - 4-digit LCD backlight display



- Memory download to PC
- Detection range: 0.00 - 0.50%
- Low battery indicator/Auto power-off

- **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.