



## Alphasport Scoop Stretcher - Metal and Aluminium

**\$783.67**

In stock  
Only **1000** left  
**SKU 99041**

 Add to Wish List  
 Add to Compare

### Details

- The Alphasport Scoop Stretcher provides superior comfort and safety to patients who need spinal immobilisation.
  - The stretcher can be separated into two parts.
  - Telescoping tubes allow the stretcher's length to be adjusted so that it can easily fit patients of various heights.
  - The open-centre design allows patients to be X-Rayed.
  - The Narrow foot end allows use in confined space immobilisation.
  - Folds compact for storage.

#### Dimensions:

Folded: 1200 x 430 x 84mm

Open: 2025 x 430 x 70mm

Weight: 8kgs

Load Limit: 250kgs

#### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all