



Blood Pressure Monitor Set Erka Switch Adult with Aneroid, Sphyg & Stethoscope

\$365.00

In stock
Only **1000** left
SKU 400142

 [Add to Wish List](#)

 [Add to Compare](#)



Details

• Features:

- Comfortable handling thanks to its flat construction
- Outstanding ergonomics and simple adaptation for left or right handed use
- Replaceable manometer housing
- Brass valve screw works smoothly and is adjustable to exact settings
- User-friendly valve cap for ease of use, even with gloves
- Latex free bladder

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

[Check items to add to the cart or select all](#)

