



Drug Safe 4 Shelves 600h x
450w x 175mm d Weight
70kg

\$1,339.67

In stock
Only **1000** left
SKU 99819

 Add to Wish List

 Add to Compare



Details

- Drug Safe 4 Shelves.
- Constructed from 10mm steel plate to Health Department specifications
- Finished in durable grey enamel paint
- Six lever pick resistant key lock
- Four 13mm mounting holdes are provided in the back plate of the safe
- Removable shelves and concealed hinges

Specifications:

- External dimensions: 450mm w x 175mm d x 600mm h
- Weight: 70kg



Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.