



Elastoplast Sensitive Dressing Length 6cm x 10cm (10)

\$15.28

In stock
Only **1000** left
SKU 2666

 Add to Wish List
 Add to Compare

Details

Elastoplast Sensitive dressing is used to protect cuts and wounds.

- Strip size is 6cm x 10cm and can be cut to desired length
- Non stick pad protects and cushions wound
- Flexible material which bends and flexes as you move
- Latex free

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all