



Fluorescein Strips (Box 100)

\$102.07

In stock
Only **994** left
SKU 99078

 [Add to Wish List](#)

 [Add to Compare](#)

Details

Sterilized applicators impregnated with fluorescein for use in ophthalmic tests. It is a fluorescent chemical (orange-red compound) that exhibits intense fluorescence in alkaline solution and is used in ophthalmology to reveal corneal lesions.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all