



Gas Air Horn Refill 250ml

\$21.52

In stock
Only **1000** left
SKU 99875

 Add to Wish List

 Add to Compare





Details

Refill for Gas Air Horn:

- Hand held and operated air horn ideal for a variety of sporting and industrial applications.
- Environmentally friendly canister does not contain any harmful gas.
- Refill 250ml

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.