



## Heine Mini 3000 2.5V Otoscope (Auroscope) Non- Fibre Optic in hard case

**\$265.53**

In stock  
Only **1000** left  
**SKU** 400069

 [Add to Wish List](#)

 [Add to Compare](#)



## Details

Heine Mini 3000 2.5V Otoscope non-fibre optic in hard case.

- Black with chrome finish
- Battery handle
- Includes five of both 2.5mm and 4mm disposable tips.

## Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

## Related Products

Check items to add to the cart or select all