



Jobst Travel Socks Beige/Black

\$46.21

In stock
Only **1000** left
SKU 330006

 Add to Wish List

 Add to Compare



Details

- Jobst Travel Socks are high-quality knee length socks, specifically designed to counteract swelling and discomfort during prolonged inactive periods of travel. The clinically proven gradient compression helps improve blood circulation in the legs, which assists in preventing swelling and aching in the legs, ankles and feet.

Features and Benefits:

- Clinically proven gradient compression
- Helps protect legs from swelling and discomfort thanks to gradient compression (15-20mmHg)
- Promotes venous blood flow in your legs
- Soft cotton blend
- Excellent breathability, leaving legs feeling fresh and comfortable
- Durable with a reinforced toe and heel
- Stylish dress sock appearance
- Ideal for men and women
- Ribbed pattern complements both business and casual attire
- Available in two fashionable colours - Beige and Black



Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all