



Mini Ice Cells Medi-Ice Pak- 32 per Sheet Single Use

€8.18

In stock
Only **949** left
SKU 310021

 Add to Wish List

 Add to Compare

Brand: Medi-Ice

- Each cell 4 x 6cm
- Non-sterile
- Disposable

Details

Medi-Ice Pak Mini Ice Cells

- 32 per sheet (4 x 6cm per cell)
- For single use only
- Perfect for fingers and toes

Store Medi-Ice in the freezer or refrigerator. Wrap your cold Medi-Ice in a light towel or cloth and apply to the injured area, moulding the pack to ensure maximum skin contact. Cold therapy can be applied and removed in 20-minute rotations for up to 48-72 hrs if merited.

The packs are easy to use and are used by health professionals, sports professionals and in first aid kits at work and



home around Australia and overseas.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all