



Qwik Hot Instant Pack 20 x 12cm Single Use

€4.88

In stock
Only **1000** left
SKU 99450

 Add to Wish List

 Add to Compare



Details

The Qwik-Hot Instant Hot Compress promotes healing by increasing oxygenation while helping to relax muscles and decrease joint stiffness. It is convenient and cost-effective. Instant Hot Compresses provide quick, localized therapy in a single-use, disposable compress.

Ingredients:

- Magnesium Sulfate
- Water

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical



emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all