



Respirator Mask N95 (50) Proshield Fluid Resistant - Made in Japan

€90.28

Out of stock

SKU 99730

 Add to Wish List

 Add to Compare

The ProShield N-95 fluid-resistant mask provides a high level of respiratory protection for clinicians/patients.

ARTG 305431

Proshield N-95 respirators comply with AS/NZS 1716:2012 Class P2, Test Safe Australia report 36251.

Proshield N-95 respirators also comply with AS4381:2015 Level 3 barrier medical face mask.

NIOSH Respirator TN01-11

TC-84A-3348 N95 protection

Details

Proshield Respirator Mask N95 (50) Fluid Resistant (50/BX)

The ProShield N-95 fluid-resistant mask provides a high level of respiratory protection for clinicians/patients who are potentially exposed to tuberculosis.

Orange single pleat particulate respirator with bands.

MADE IN JAPAN

Features and Benefits:



Bacterial Filtration Efficiency (BFE)

With a Bacterial Filtration Efficiency (BFE) of greater than 99% for media greater than 3 microns, this mask provides superior protection from potential infection.

Particulate Filtration Efficiency (PFE)

The super high Particulate Filtration Efficiency (PFE) material filters >99% of particles greater than 0.1 microns. When securely fitted, the mask helps to provide an effective barrier to particulate matter from laser plumes and electrocautery procedures.

PFE measures effectiveness of the mask material to filter particulate matter.

Particulate Filter Respirator

The N-95 mask is NIOSH approved as a N-95 particulate filter respirator. It meets or exceeds the standard performance criteria demanded by NIOSH for the management of Tuberculosis (TB)

NIOSH - National Institute of Occupational Health and Safety (USA).

Fluid Resistance

The fluid resistant qualities of the mask provide protection against fluid strikethrough.

Fluid resistance is measured in accordance with ASTM F1862-98 & ASTM F1862-00a standard test method for resistance of medical face masks to penetration by synthetic blood at 16.0kPa pressure (120mmHg).

ARTG 305431

Proshield N-95 respirators comply with AS/NZS 1716:2012 Class P2, Test Safe Australia report 36251.

Proshield N-95 respirators also comply with AS4381:2015 Level 3 barrier medical face mask.

NIOSH Respirator TN01-11

TC-84A-3348 N95 protection

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



More Information

SPECIAL PRICE

€104.72

Frequently bought together

You may also need



[Facial Tissues \(224\) - Single Box](#)

€3.95

[Add to cart](#)



[Sharps Container Round 1.8L Screw Lid with Antispill Insert](#)

€10.73

[Add to cart](#)



[Versa Towel Large 49.5 x 41.5cm - Single](#)

Special Price €20.72

Regular Price €22.59

[Add to cart](#)



[Yellow Biohazard Waste Bags 720 x 58mm \(10\) 40L](#)

€6.10

[Add to cart](#)